360° Head Coverage

360° KOPFABDECKUNG

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PRODUCT DETAILS

Materials:

High elastic lycra + malleable gel

Ingredients:

30% water, 48% glycerin and 20% Potassium Polyacrylate

Applicable Head Circumferences:

42-62 cm

Product size:

17 x 22 x 16 cm

Color:

Black

Weight:

600 g

Cold or heat therapy

Improve your sleep quality.
Can relieve migraines.

Put in the freezer for 1.5 - 2 hours

Microwave for 30 seconds or more

INSTRUCTION

Compared with cold therapy, some people prefer relaxing themselves by the hot one, while our headache relief cap can apply to these two therapies at the same time.

For cold Compress: Freeze the cap in a refrigerator for 2 hours, and then take it out and put it on to relieve your migraine or facial tension. For hot compress: Heat the cap in a microwave (700W) for 40 to 60 seconds, and then take it out and put it on to relieve your migraine or facial tension.

Notice:

1. The cap may be too cold after a long time of refrigeration, it's recommended to put a towel or cushion as a buffer on your ears and eyes to avoid frostbite.
2. After heating, please apply the hot compress way according to your own suitable temperature. If you feel uncomfortable, please take it off immediately.
3. This product can NOT be washed in water, but you can wipe it with towels.
4. Due to the special material of the cap, if the cap has a slight smell when you received, please put it in a ventilation a few days to remove the smell.
5. Do NOT use the product directly on the wounded skin.

Storage:

Keep the product in a dry and ventilated place, avoiding sun exposure.

GUARANTEE/WARRANTY

The warranty period is 2 years from the date of purchase. This includes a 12-month guarantee from the date of purchase. This guarantee is valid for the whole of Europe. The purchase receipt is regarded as proof of purchase. The customer undertakes to read the instructions carefully before making a claim under the guarantee/warranty in order to rule out any operating errors.

The product has been thoroughly checked for quality and functionality before dispatch. Exclusion of the guarantee occurs in the following cases:

- in the event of changes to the product without our approval
- in the event of use other than for the intended purpose
- in the event of defects caused by improper maintenance and/or repair
- in the event of commercial use of the article
- in the event of damage caused by falling

When making a claim under the guarantee/warranty, please return the complete product including accessories in a packaging safe for transport to the point of sale from which the product was purchased. Describe the existing defect in the product and enclose the original proof of purchase with the shipment.

The seller reserves the right to subsequent performance before refunding the purchase price. Necessary transport costs will be reimbursed to the customer if a justified claim of the legal warranty/guarantee is established. The exercise of the statutory rights in respect of defects shall be free of charge. These rights are not limited by the guarantee.

SHP International Trading GmbH (Rathausstraße 5, 96342 Stockheim, Germany) hereby declares that this product is in conformity with the essential requirements of the applicable European Directives.

Intended use

This product (360° Head Coverage) is only intended for use as a Head Coverage. Any other use or use that goes beyond this is deemed to be improper.

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